

Column: **MEGA-SE54 FAST - 0.10mm, 0.10 μ m, 5m**

Catalog Code: F-SE54-010-010-5

Conditions:

Injection: Split 230°C, 0.5 μ L, 1:300 Split Ratio.

Sample dilution: 1% in Cyclohexane.


Detector: FID, 250°C.

Oven Temp.: 50°C (0.1 min), 50°C/min, 250°C (5 min).

Carrier Gas: Hydrogen, 0.5 mL/min.

Peaks:

- | | |
|------------------------|--|
| 1. Benzyl alcohol | 17. Alpha-isomethyl ionone |
| 2. Limonene | 18. Alpha-methyl ionone |
| 3. Linalool | 19. Lilial |
| 4. Veratrol | 20. / |
| 5. Methyl 2-octanoate | 21. Lyral 1 + Lyral 2 |
| 6. Citronellol | 22. Farnesol 1 |
| 7. Citral 1 | 23. Farnesol 2 |
| 8. Geraniol | 24. Farnesol 3 |
| 9. Cinnamic aldehyde | 25. Amyl Cinnamal |
| 10. Citral 2 | 26. Hexyl cinnamal + Benzyl Salicylate |
| 11. Anisel Alcohol | 27. 1-Phenyl-decanone |
| 12. Hydroxycitronellal | 28. / |
| 13. Cinnamic Alcohol | 29. Benzyl cinnamate |
| 14. Eugenol | |
| 15. Coumarine | |
| 16. isoeugenol | |



Acknowledgement: Prof C. Bicchi, C. Brunelli et al. - Università di Torino, Dipartimento di Scienza e Tecnologia del Farmaco, Via Pietro Giuria, 9 - Torino - Italy.