

2-phenylpropanol - enantiomers


Chiral separation on MEGA-DEX columns

CHIRAL DEX

MEGA[®]
gas chromatography solutions

Column: **MEGA-DEX B-SE - 0.25mm, 0.25 μ m, 30m**

Catalog Code: DEX-BSE-025-025-30



Acknowledgments: author of this work is Dr. A. Adler, University of Köln, Organic Chemistry Dept., Greinstrasse 4, Köln, Germany. Many thanks also to and Dr. R. Forche, "Dr. R. Forche Chromatographie" (www.gc-hplc.de).