

Acetic Acid monitoring in Ethyl Acetate



Column: **MEGA-FFAP EXT FAST - 0.15mm, 0.25 μ m, 10m**

Catalog Code: F-FFAPEXT-015-025-10

Conditions:

Injection: Split 250°C, 1.5 μ L, 25mL/min split flow.

Oven Temp.: 35°C (0.5min), 20°C/min.

Detector: FID, 250°C.

Carrier Gas: Helium, 100kPa.

• = Ethanol impurities and other residual solvents

