

Biodiesel - Glycerin mono-, di- and Triglycerides (UNI EN ISO 14105, ASTM 6584)

CUSTOM
DEDICATED

MEGA®
gas chromatography solutions

Column: **MEGA-BIODIESEL 105 - 0.32mm, 0.10µm, 15m**

Catalog Code: C-BIODIESEL105-032-010-15

Peaks:

1. Glycerin
2. I.S. 1
3. Monoolein
4. I.S. 2
5. Diolein
6. Triolein

Conditions:

Injection: PTV Program Temp. 50°C, 999°C/min, 230°C (15min).

Sample: Standards derivatized by MSTFA.

Oven Temp: 50°C (1min), 15/min, 180°C, 7°C/min, 230°C, 10°C/min, 370°C (5min).

Detector: FID 380°C.

Carrier Gas: Helium, 3 mL/min.

