

Champor

Chiral separation on MEGA-DEX columns

CHIRAL DEX

MEGA[®]
gas chromatography solutions

Phase: **MEGA-DEX DET Beta**
Catalog code: DEX-DET-B-025-025-25

Phase: **MEGA-DEX DMT Beta**
Catalog code: DEX-DMT-B-025-025-25

Phase: **MEGA-DEX DAC Beta**
Catalog code: DEX-DAC-B-025-025-25


Columns dimensions: 0.25mm, 0.25 μ m, 25m.

Conditions:

Injection: Split 250°C, 1.5 μ L, 1:50 Split Ratio.

Sample dilution: 1mg/1mL in n-Hexane.

Contact us for additional information.



Acknowledgement: Prof C. Bicchi, C. Cordero et al. - Università di Torino, Dipartimento di Scienza e Tecnologia del Farmaco, Via Pietro Giuria, 9 - Torino - Italy.