

FAMES C6 - C24 (Fatty Acid Methyl Esters)

STANDARD



Peaks:

1. C 6:0	11. C 16:0	21. C 18:4 n-3	31. C 22:2 n-6
2. C 8:0	12. C 16:1 n-7	22. C 20:0	32. C 24:0
3. C 10:0	13. C 17:0	23. C 20:1 n-9	33. C 22:6 n-3
4. C 11:0	14. C 17:1 n-7	24. C 20:2 n-6	34. C 24:1 n-9
5. C 12:0	15. C 18:0	25. C 20:3 n-6	
6. C 13:0	16. C 18:1 <i>cis</i> n-7	26. C 20:4 n-6	
7. C 14:0	17. C 18:1 <i>trans</i> n-7	27. C 20:3 n-3	
8. C 14:1 n-5	18. C 18:2 n-6	28. C 20:5	
9. C 15:0	19. C 18:3 n-6	29. C 22:0	
10. C 15:1 n-5	20. C 18:3 n-3	30. C 22:1 n-9	

Column: **MEGA-WAX - 0.32mm, 0.25 μ m, 25m**

Catalog Code: S-WAX-032-025-25

Conditions:

Injection: On-Column.

Oven Temp: 60°C (2min), 15°C/min, 105°C, 3°C/min, 195°C, 5°C/min, 220°C (1min).

Detector: FID, 250°C.

Carrier Gas: Hydrogen, 30 kPa.

