

Aqueous residual from phenolic resins production plant



Column: **MEGA-WAX Plus - 0.25mm, 0.25 μ m, 30m**

Catalog Code: C-WAXP-025-025-30

Conditions:

Injection: Split, 250°C.

Detector: FID, 250°C.

Oven Program: 150°C, 8°C/min.

Carrier Gas: Hydrogen, 60kPa.

Peaks:

1. Methanol (used as solvent)
2. unknown
3. Dimethylphenol (Internal Standard, 1% in Methanol)
4. Hexamine (1% in Methanol)
5. Phenol

The rest of the sample is aqueous.

