

Benzene and aromatics volatile compounds in air

STANDARD

MEGA[®]
gas chromatography solutions

Column: **MEGA-WAX - 0.25mm, 0.25 μ m, 30m**

Catalog Code: S-WAX-025-025-30

Conditions:

Injection: Thermal Desorber (Tenax GR trap 300°C),
Split Injector 80°C, 10mL/min.

Detector: FID, 200°C.

Oven Program: 40°C (5min), 5°C/min, 200°C.

Carrier Gas: Helium, 1.5 mL/min, constant flow.

Peaks:

1. Benzene
2. Toluene
3. Ethylbenzene
4. p-Xylene
5. m-Xylene
6. o-Xylene

