


# Isoborneol (+/-)

Chiral separation on MEGA-DEX columns

CHIRAL DEX

MEGA<sup>®</sup>  
gas chromatography solutions



Phase: **MEGA-DEX DET Beta**  
Catalog code: DEX-DET-B-025-025-25

Phase: **MEGA-DEX DMT Beta**  
Catalog code: DEX-DMT-B-025-025-25

Phase: **MEGA-DEX DAC Beta**  
Catalog code: DEX-DAC-B-025-025-25

Phase: **MEGA-DEX DMP Beta**  
Catalog code: DEX-DMP-B-025-025-25

Columns dimensions: 0.25mm, 0.25  $\mu$ m, 25m.

Acknowledgement: Prof C. Bicchi, C. Cordero et al. - Università di Torino, Dipartimento di Scienza e Tecnologia del Farmaco, Via Pietro Giuria, 9 - Torino - Italy.