

Purgeable Aromatics (EPA Method 602)

STANDARD

MEGA[®]
gas chromatography solutions

Column: **MEGA-WAX - 0.53mm, 1.00 μ m, 25m**

Catalog Code: S-WAX-053-100-25

Conditions:

Injection: Split 250°C, 1.5 μ L.

Sample dilution: 2000 μ g/mL each in P&T Methanol.

Detector: FID, 250°C.

Oven Temp.: 50°C, 8°C/min, 120°C.

Carrier Gas: Helium, 20 kPa.

Peaks:

1. Benzene
2. Toluene
3. Ethylbenzene
4. Chlorobenzene
5. 1,3-dichlorobenzene
6. 1,4-dichlorobenzene
7. 1,2-dichlorobenzene

