

Nitrogen/Phosphorus Pesticides & Herbicides - Triazines (EPA method 619)

STANDARD

MEGA[®]
gas chromatography solutions

Column: **MEGA-WAX - 0.32mm, 0.25 μ m, 25m**

Catalog Code: S-WAX-032-025-25

Conditions:

Injection: Split 250°C, 1.5 μ L, 1:70 Split Ratio.

Detector: FID, 300°C.

Oven Temp.: 230°C isothermal.

Carrier Gas: Helium, 100 kPa.

Peaks:

Peaks:	CAS #:
1. Trietazine	1912-26-1
2. Prometon	1610-18-0
3. Propazine	139-40-2
4. Terbutylazine	5915-41-3
5. Atrazine	1912-24-9
6. Prometryne	7287-19-6
7. Terbutryne	886-50-0
8. Simazine	122-34-9
9. Ametryne	834-12-8
10. Simetryne	1014-70-6

