


Bergamot Oil



Column: **MEGA-1701 FAST - 0.10mm, 0.10 μ m, 5m**

Catalog Code: F-1701-010-010-5



Conditions:

Injection: Split 230°C, 0.5 μ L, 1:250 Split Ratio.

Sample dilution: 1% in Cyclohexane.

Oven Temp: 50°C (0.1min), 50°C/min, 250°C.

Detector: FID 250°C

Carrier Gas: Hydrogen, 0.5 mL / min.

Peaks:

1. Alpha-Pinene
2. Beta-Pinene
3. Myrcene
4. Limonene
5. p-Cymene
6. Gamma-Terpinene
7. Linalool
8. Linalyl Acetate

Acknowledgments: Prof C. Bicchi, C. Brunelli et al., Università di Torino, Dipartimento di Scienza e Tecnologia del Farmaco, Via Pietro Giuria, 9 Torino, Italy.