


Dimethylanilines

STANDARD

MEGA[®]
gas chromatography solutions



Column: **MEGA-WAX - 0.32mm, 0.25 μ m, 30m**

Catalog Code: S-WAX-032-025-30

Conditions:

Injection: Split 250°C, 1 μ L.

Oven Temp: 130°C isothermal.

Detector: FID 250°C.

Carrier Gas: Nitrogen, 1.8 mL/min.

Peaks:

1. 2,6 dimethylaniline
2. 2,4 dimethylaniline
3. 2,5 dimethylaniline
4. 3,5 dimethylaniline
5. 2,3 dimethylaniline
6. 3,4 dimethylaniline