

# Alkanes C<sub>1</sub> - C<sub>4</sub>, CO<sub>2</sub> and H<sub>2</sub>S

CUSTOM  
DEDICATED

MEGA®  
gas chromatography solutions

Column: **MEGA-PS264 - 0.32mm, 5.00μm, 60m**

Catalog Code: C-PS264-032-500-60

**Conditions:**

Injection: Split 180°C, 100μl injection volume, 50 mL/min Split flow.

Detector: TCD, 180°C.

Oven Temp.: 40°C isothermal.

Carrier Gas: Hydrogen, 65 kPa.

**Peaks:**

1. CH<sub>4</sub>
2. CO<sub>2</sub>
3. Ethane
4. H<sub>2</sub>S
5. N-propane
6. Isobutane
7. N-butane

